



HealkerenTM

DIAMOND TIPS



Diamond tip comes usually with wands of different sizes to provide an ideal resurfacing depth and to better suit different parts of the face and body. Natural diamond tips are attached to the tip of the wand to remove the dead cells and they are vacuumed away.

M/s. Healkeren Medtronics Private Limited

No.4, 2nd Cross, Lakshmipura Main Road, Krishnappa Industrial Estate, Abbigere, Bengaluru - 560 090, Karnataka, India.

Mobile : 93101 18774 / 80882 84888 / 80882 84898, Email : sales@healkeren.com / healkeren@gmail.com / info@healkeren.com

www.healkeren.com