



Now Real 'PICO' Laser!

# Pico Premium





**01 Stable Operating**

- Minimal pain

**02 Convenient to Treat**

- Smart protocols

**03 Low Risk of Side Effect**

- Fast treatment

■ Advantage

- Stable Operation
- Convenient to Treat with Smart Protocols
- Low risk of Side Effect
- Fast Treatment

■ The Real PICO

Thanks to higher peak power, PICO can treat all kinds of tattoos and pigment lesions more faster and easier.

■ Pico Premium has the shortest pulse width of 350 ps

■ Low risk of side effect

- Less thermal effect on the treated area makes the best result of treatment with low risk of side effect


■ Indications

Tattoo, Wrinkle, Acne, Nevus of ota, Blue nevus, Speckled nevus, Lentigo, Cafe au lait, Hyperpigmentation, Freckle, Melasma and so on.

■ Super ultra short pulsed width

Picosecond Nd:YAG Laser (1064nm & 532nm + 595nm 660nm HP+ Fractional HP) with super ultra short pulse width

■ Hand Piece



Fractional Hand Piece for 1064nm (Optional)

Zoom Hand Piece (2-10mm, Standard)

Dye Hand Piece (585nm, Optional)

Dye Hand Piece (650nm, Optional)

Dye Hand Piece (595nm, Optional)

Dye Hand Piece (660nm, Optional)



Healkeren<sup>TM</sup>

Before and After



SPECIFICATIONS

Model	PICO PREMIUM	
Wavelength	1064nm/532nm	
Operating Mode	585nm/650nm/595nm/660nm	
Pulse Energy	1064nm	Max. 500mj
	532nm	Max. 250mj
	585/595nm	Max. 150mj
	650/660nm	Max. 100mj
Pulse Width	Nd:YAG	350ps
Spot Size	2~10mm (1064/532nm)	
	2mm-3mm (55/650nm/595/660)	
Pulse Rate	2Hz-10Hz	
	585/650 : 1-2 Hz	
Aiming Beam	Diode 650nm	
Optional Handpiece	- Dye Hand Piece (585, 650nm / 595, 660nm) - Fractional Handpiece	

**M/s. Healkeren Medtronics Private Limited**

No.4, 2nd Cross, Lakshmipura Main Road, Krishnappa Industrial Estate, Abbigere, Bengaluru - 560 090, Karnataka, India.

Mobile : 93101 18774 / 80882 84888 / 80882 84898, Email : sales@healkeren.com / healkeren@gmail.com / info@healkeren.com

[www.healkeren.com](http://www.healkeren.com)