



IDS CRYO-ZET

CE

AIR COOLING SYSTEM

Cryo-Zet





The most powerful
Air Cooling System



01

Convenient and economical for users

- It is very economical because treatment times can be considerably reduced.
- It is very convenient because using is very easy.
- It is very useful because consumable parts are not needed.

02

Very effective and safe for patients

- It can reduce skin damage because it can be used with lasers.
- It can reduce pain and recovery time.
- It can protect skin because cooling can be made evenly.

03

The suitable lasers

- Fractional CO2 Laser, IPL, RF Microneedle type, Er YAG, Diode laser, Alexandrite, Nd:YAG(Q-Switched, Long Pulsed YAG)

04

Special Features

- ✓ Air flow controlling system.
- ✓ Cooling air temperature (Max. -35 degree).
- ✓ Quick cooling with high quality compressor.
- ✓ Auto adjustable controlling system.
- ✓ Auto detecting system to keep best condition.

SPECIFICATIONS

SPECIFICATIONS	
A.C Input	220~240V, 50/60Hz
Air Level	10 Steps
Size	400(W)x610(L)x963(H)
Air Flow Output	1000L/Min
Temperature	0~-35°C
Weight	50kg

M/s. Healkeren Medtronics Private Limited

No.4, 2nd Cross, Lakshmipura Main Road, Krishnappa Industrial Estate, Abbigere, Bengaluru - 560 090, Karnataka, India.

Mobile : 93101 18774 / 80882 84888 / 80882 84898, Email : sales@healkeren.com / healkeren@gmail.com / info@healkeren.com

www.healkeren.com