

OXYGEN HOOD : Heal - OXN/OSP



ISO 9001:2008
ISO 13485:2003



Healikeren™

OXYGEN HOOD : Heal - OXN/EXP

Heal-OXN and Heal-EXP oxygen hoods are transparent polycarbonate shell that surrounds the baby head. These hoods can be used effectively to avoid direct flow of cold oxygen on baby's face and head. The shell is designed in such a way so that other body parts like chest, stomach etc., can be accessed. The opening is generally softly padded keeping in mind the soft skin of the baby so that they don't hurt themselves. The oxygen hood is available in two sizes.

SILENT FEATURES

Heal-OXN Size Height 6.5" x Width 8.5"

- Single piece round shape made up of Autoclavable polycarbonate material
- Trauma free silicone neck adjustment flap
- Bilateral oxygen nozzle prevents direct flow of cold oxygen on patient's head
- 2 meter oxygen tubing

Heal-EXP Size Height 7.5" x Width 8.5"

- Single piece round shape made up of Autoclavable polycarbonate material
- Trauma free silicone neck adjustment flap
- Bilateral oxygen nozzle prevents direct flow of cold oxygen on patient's head
- 2 meter oxygen tubing

M/s. Healkeren Medtronics Private Limited

No.4, 2nd Cross, Lakshmipura Main Road, Krishnappa Industrial Estate, Abbigere, Bengaluru - 560 090, Karnataka, India.

Mobile : 93101 18774 / 80882 84888 / 80882 84898, Email : sales@healkeren.com / healkeren@gmail.com / info@healkeren.com

www.healkeren.com